

Before & After



Courtesy of G. Mercik, MD



Courtesy of D. Day, MD



Courtesy of D. Day, MD

Proven results and highly satisfied patients

Several clinical studies of HAIRestart® have demonstrated high efficacy. Significant improvement of hair growth along with increased hair density and a stabilization of hair loss were reported. Patients typically report **more elastic and less dry scalp**, and **thicker, more durable, healthier and denser hair** than before.

HAIRestart®

Non-invasive laser stimulation of hair growth



Fotona
choose perfection



What is HAIRestart®?

Fotona's HAIRestart® laser hair-growth solution is an effective method of treating hair loss disorders in men and women. This innovative solution is based on gentle, non-invasive laser stimulation of dormant hair follicles to **regenerate the scalp and promote hair growth**.

How does it work?

HAIRestart® takes advantage of the unique, patented Fotona SMOOTH® laser modality, delivering precisely engineered **laser pulses** to the tissue, triggering the scalp in a safe and therapeutic way with ultimate comfort for the patient. The result is **increased tissue regeneration** and the promotion of **new follicle growth**.



More About HAIRestart®

Key advantages of HAIRestart®

- Prevents further hair loss
- Increases hair density and thickness
- Promotes new follicle growth
- No medications or down-time

What is the treatment like?

The HAIRestart® procedure is **well tolerated by patients**, with minimal discomfort. There are several passes of the laser over the scalp during the treatment and patients feel only a slow and gradual heat build-up. There is some redness, which subsides in a few days, and there can be some slight delayed peeling of the skin. The duration of the treatment depends on the size of the affected area. In general it is about a 30 minute appointment.

Before & After

